ZAA 2300 Automatic film applicator coater

- ▶ Automatic laboratory equipment for accurate and reproducible application of coating materials, adhesives and similar products, independent of the operator
- ▶ Multifunctional use with reversible, double sided glass plate (glass/printing blanket), easy to turn over or to change
- ▶ Adjustable application area with moveable start and stop positions
- ▶ For use with almost all types of film applicators
- ▶ Optional vacuumplate
- ▶ Optional versions equipped with isolation plate for use with heatable vacuum plates or heating plates
- ▶ Easy to handle
- ▶ Reliable results
- ▶ Modern and functional design

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Manual applications may have irregularities caused by uneven pressure or non-uniform drawing speed. The ZAA 2300 Automatic film applicator coater ensures an accurate and reproducible coating.

Features

- Double sided glass plate for dual-purpose operation: one side with printing blanket for wire-bar applicators and profile rods, on the other side glass surface for other applications
- Option: vacuumplate for substrates; the vacuum created fixes the substrate evenly in place

Standard delivery

- 1 automatic film applicator coater
- 1 reversible, double sided glass plate / printing blanket
- 1 polyester mat
- 1 spirit level
- 1 certificate of manufacturer
- 1 instruction manual

Options

- ACC121 precision vacuumplate
- ACC039/ACC122 (230V/115V) vacuum pump with hose
- ACC225.230/ACC225.115 (230V/115V) heating plate /
- ACC188/ACC189 (230V/115V) heated precision vacuum plate, all for temperature up to 150°C (302°F), those options can only be used together with an automatic film applicator coater equipped with an isolation plate (2300H/FH/FFH)
- ZTC 2200 Test charts
- Film applicators with maximum outer width of 300 mm (11.81")

Handling

- Choose the applicator(s).
- Choose the appropriate base (glass surface, printing blanket for wire bar applicators and profile rods or vacuumplate).
- Mount wire bar applicator or profile rod to the ZAA 2300 Automatic film applicator coater and fix with the weight support.
- Position other applicators in front of the weight support.
- Adjust the desired application length.
- Fix the desired substrate, e.g. test chart, with the start position clamp.
- Set the desired drawing speed.
- Pour the product to be tested in front of the applicator in drawing direction.
- Start the ZAA 2300 Automatic film applicator coater.
- After the application remove the applicator, the weight and the substrate.

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